



Product Datasheet

PYCR2 antibody (orb667861)


Description	Rabbit polyclonal antibody to PYCR2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human PYCR2. The exact sequence is proprietary.
Target	PYCR2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q6AY23 , Q922Q4 , Q96C36
Entrez	69051 , 29920 , 364064
Dilution Range	WB: 1-500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of PYCR2 expression...



Immunohistochemical analysis of PYCR2 staining...



Immunofluorescent analysis of PYCR2 staining...