

Product Datasheet

IL33 antibody (orb66787)

Description	Rabbit polyclonal to IL33. Interleukin 33 (IL33) is a cytokine that binds to and signals through IL1R...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	N-terminal region of human IL33
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	31
Uniprot ID	Q8BVZ5
Dilution Range	ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.