


Product Datasheet

DNER antibody (orb66788)

Description	Rabbit polyclonal to DNER. Delta/notch-like EGF repeat containing (DNER) is a single pass type I memb...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human DNER
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	78
Uniprot ID	Q8NFT8
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...