



Product Datasheet

Ribosomal Protein L21 antibody (orb667985)


Description	Rabbit polyclonal antibody to Ribosomal Protein L21
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human Ribosomal Protein L21 . The exact sequence is proprietary.
Target	RPL21
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P20280, O09167, P46778
Entrez	6144, 79449
Dilution Range	WB: 1-500-2000, IHC-P: 1-50-200, IF: 1-50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Ribosomal Prote...



Immunohistochemical analysis of Ribosoma...



Immunofluorescent analysis of Ribosomal

...