


Product Datasheet

KLF4 antibody (orb66809)

Description	Rabbit polyclonal to KLF4. Kruppel-like factor 4 (gut) (KLF4) is a transcription factor which acts as...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human KLF4
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	55
Uniprot ID	O43474
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...