



Product Datasheet

SLC18A3 antibody (orb668126)

Description	Rabbit polyclonal antibody to SLC18A3
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human SLC18A3. The exact sequence is proprietary.
Target	SLC18A3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q16572 , Q62666 , O35304
Entrez	6572 , 60422
Dilution Range	WB: 1-500-2000
Expiration Date	12 months from date of receipt.




Western blot analysis of SLC18A3 express...



Immunohistochemical analysis of SLC18A3

...



Immunofluorescent analysis of SLC18A3 st...