

Product Datasheet

MMP13 antibody (orb66815)

Description	Rabbit polyclonal to MMP13.Matrix metalloproteinase 13 (collagenase 3) (MMP13) is a zinc-dependent en...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	C-terminal region of human MMP13
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	54
Uniprot ID	P45452
Dilution Range	ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.