

Product Datasheet

CXCL16 antibody (orb66816)

Description	Rabbit polyclonal to CXCL16.CXCL16 acts as a scavenger receptor on macrophages which specifically bi...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	N-terminal region of mouse CXCL16
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	28
Uniprot ID	Q8BSU2
Dilution Range	ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.