


Product Datasheet

GMCSF antibody (orb66817)

Description	Rabbit polyclonal to GMCSF. Granulocyte-macrophage colony-stimulating factor (GMCSF) is also known as...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	N-terminal region of mouse GMCSF
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	16
Uniprot ID	P0158
Dilution Range	WB:1:1000-1:5000, ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant mouse GMCSF.