


Product Datasheet

CCL4 antibody (orb66821)

Description	Rabbit polyclonal to CCL4. Chemokine (C-C motif) ligand 4 (CCL4) is a monokine with inflammatory and ...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	C-terminal region of mouse CCL4
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	10
Uniprot ID	P14097
Dilution Range	WB:1:1000-1:5000, ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant mouse CCL4.