


Product Datasheet

CCL22 antibody (orb66826)

Description	Rabbit polyclonal to CCL22.Chemokine (C-C motif) ligand 22 (CCL22) is a cytokine with chemotactic ac...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	N-terminal region of mouse CCL22
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	11
Uniprot ID	O88430
Dilution Range	WB:1:1000-1:5000,ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant mou...