


Product Datasheet

FGF2 antibody (orb66832)

Description	Rabbit polyclonal to FGF2. Fibroblast growth factor 2 (basic) (FGF2) is a member of the FGF family of...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	N-terminal region of mouse FGF2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	31
Uniprot ID	Q7TPG9
Dilution Range	WB:1:1000-1:5000, ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...