

Product Datasheet

Tachykinin Receptor 2 antibody [Out of stock] (orb668336)

Description	Rabbit polyclonal antibody to Tachykinin Receptor 2
Species/Host	Rabbit
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Tachykinin Receptor 2. The exact sequence is proprietary.
Target	TACR2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P21452
Entrez	6865
Dilution Range	WB: 1-500-2000
Expiration Date	12 months from date of receipt.