


## Product Datasheet

### EHMT1 antibody (orb66842)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to EHMT1.EHMT1 is a Histone-lysine N-methyltransferase that specifically mono-and ... |
| <b>Species/Host</b>        | Rabbit                                                                                                  |
| <b>Reactivity</b>          | Human                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                            |
| <b>Tested Applications</b> | WB                                                                                                      |
| <b>Immunogen</b>           | C-terminal region of human EHMT1                                                                        |
| <b>Preservatives</b>       | 0.01% NaN <sub>3</sub>                                                                                  |
| <b>Storage</b>             | -20°C                                                                                                   |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                              |
| <b>Purity</b>              | Affinity purification                                                                                   |
| <b>MW</b>                  | 141                                                                                                     |
| <b>Uniprot ID</b>          | <a href="#">Q9H9B1</a>                                                                                  |
| <b>Dilution Range</b>      | WB:1:1000-1:5000                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                         |



Western blot analysis of total cell extr...