


Product Datasheet

RPS6KA3 antibody (orb66843)

Description	Rabbit polyclonal to RPS6KA3. Ribosomal protein S6 kinase, 90kDa, polypeptide 3 (RPS6KA3) is a member o...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human RPS6KA3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	84
Uniprot ID	P51812
Dilution Range	WB: 1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...