



Product Datasheet

TMED3 antibody (orb668433)


Description	Rabbit polyclonal antibody to TMED3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human TMED3. The exact sequence is proprietary.
Target	TMED3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q78IS1 , Q9Y3Q3 , Q6AY25
Entrez	300888 , 66111 , 23423
Dilution Range	WB: 1-500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of TMED3 expressio...



Immunohistochemical analysis of TMED3 st...



Immunofluorescent analysis of TMED3 stai...