


Product Datasheet

ING1 antibody (orb66846)

Description	Rabbit polyclonal to ING1. Inhibitor of growth protein 1 (ING1) cooperates with p53/TP53 in the negat...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human ING1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	47
Uniprot ID	Q9UK53
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell extr...