


Product Datasheet

AURKB antibody (orb66849)

Description	Rabbit polyclonal to AURKB. AURKB is a key regulator for the onset of cytokinesis during mitosis and ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human AURKB
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	39
Uniprot ID	Q96GD4
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell extr...