


Product Datasheet

MYST1 antibody (orb66857)

Description	Rabbit polyclonal to MYST1.MYST histone acetyltransferase 1 (MYST1) belongs to the MYST family of hi...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MYST1
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	52
Uniprot ID	Q9H7Z6
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...