

Product Datasheet

Pectinesterase antibody (orb6686)

Description	Pectinesterase antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from Arabidopsis thaliana (thale cress) pectinesterase inhibitor 18 Pectinesterase (451-557/557 aa)
Target	Pectinesterase inhibitor 18
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	Q1JPL7
Expiration Date	12 months from date of receipt.