


Product Datasheet

PRMT3 antibody (orb66860)

Description	Rabbit polyclonal to PRMT3. Protein arginine methyltransferase 3 (PRMT3) is a type I protein arginine...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal propeptide region of human PRMT3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	48
Uniprot ID	O60678
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...