


Product Datasheet

IL17F antibody pair (orb66865)

Description	Rabbit polyclonal to IL17F. 50 μl IL17F (N-term) &50 μl IL17F (C-term1). vPair antibodies repr...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	-
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	18
Uniprot ID	Q96PD4
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant hum...