

Product Datasheet

IL21 antibody (orb66866)

Description	Rabbit polyclonal antibody to IL21
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	C-terminal region of human IL21
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	18
Uniprot ID	Q9HBE4
Dilution Range	ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.