


Product Datasheet

CD86 antibody (orb66871)

Description	Rabbit polyclonal to CD86.T-lymphocyte activation antigen CD86 is a receptor involved in the costimu...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human CD86
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	38
Uniprot ID	P42081
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...