

Product Datasheet

MFGE8 antibody (orb66872)

Description	Rabbit polyclonal to MFGE8. Milk fat globule-EGF factor 8 protein (MFGE8) plays a role in the mainten...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MFGE8
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	48
Uniprot ID	Q08431
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.

