


Product Datasheet

SMYD3 antibody (orb66873)

Description	Rabbit polyclonal to SMYD3.SET and MYND domain containing 3 (SMYD3) is a histone methyltransferase t...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human SMYD3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	49
Uniprot ID	Q9H7B4
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...