


Product Datasheet

NCOA1 antibody (orb66875)

Description	Rabbit polyclonal to NCOA1. NCOA1 is a nuclear receptor coactivator that directly binds nuclear recep...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human NCOA1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	157
Uniprot ID	Q15788
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...