

Product Datasheet

SPNS2 antibody (orb668768)

Description	SPNS2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human SPNS2. Please contact us for the exact immunogen sequence. The peptide is available as orb668773.
Target	SPNS2
Form/Appearance	10 mM PBS, 0.02% sodium azide, 10 mg/ml rAlbumin
Storage	Store at 4°C.for up to two weeks. For long term storage, aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	Q8IVW8
Dilution Range	IHC-P: 1:50~500, WB:1:100~1000
Expiration Date	12 months from date of receipt.