

Product Datasheet

CIC-7 antibody (orb668772)

Description	CIC-7 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CIC-7. Please contact us for the exact immunogen sequence. The peptide is available as orb668777.
Target	CIC-7
Form/Appearance	10 mM PBS, 0.02% sodium azide, 10 mg/ml rAlbumin
Storage	Store at 4°C.for up to two weeks. For long term storage, aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	P51798
Expiration Date	12 months from date of receipt.