

Product Datasheet

PD-L1 blocking peptide (orb668775)

Description	PD-L1 blocking peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	PD-L1
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb668770 (PD-L1 antibody)
Expiration Date	6 months from date of receipt.