


Product Datasheet

KITL antibody (orb66882)

Description	Rabbit polyclonal to KITL. KIT ligand (KITL) is the binding partner of the tyrosine-kinase receptor K...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	C-terminal region of human KITL
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	31
Uniprot ID	P21583
Dilution Range	WB:1:1000-1:5000, ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...