

Product Datasheet

IL5 antibody (orb66883)

Description	Rabbit polyclonal to IL5. Interleukin 5 (IL5) is a cytokine that acts as a growth and differentiation...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	N-terminal region of mouse IL5
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	15
Uniprot ID	P04401
Dilution Range	ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.