


Product Datasheet

KRT19 antibody (orb66885)

Description	Rabbit polyclonal to KRT19. Keratin 19 (KRT19) is a member of the keratin family. Keratins are interm...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human KRT19
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	44
Uniprot ID	P08727
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...