

Product Datasheet

KRT19 antibody pair (orb66886)

Description	Rabbit polyclonal to KRT19. 50 μl KRT19 (C-term)&50 μl KRT19 (N-term).vPair antibodies represe...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	-
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	44
Uniprot ID	P08727
Expiration Date	12 months from date of receipt.

