

Product Datasheet

Human CD59 Protein (orb668873)

Description	CD59, also known as membrane attack complex inhibition factor (MACIF), Protectin, 1F5 antigen, HRF-2...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD59
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 10.8 kDa. The protein migrates as 15-18 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Leu 26 - Asn 102
Purity	95%
Source	Human CD59, His Tag (orb668873) is expressed from human 293 cells (HEK293). It contains AA Leu 26 - Asn 102 (Accession # P13987-1).
MW	10.8 kDa
Expiration Date	6 months from date of receipt.