

Product Datasheet

Rabbit LDL R Protein (orb668886)

Description	Low-Density Lipoprotein (LDL) Receptor is also known as LDLR, FH, FHC, LDLCQ2, and is a mosaic prote...
Reactivity	Rabbit
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	LDL R
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 109.2 kDa. The protein migrates as 120-140 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Ala 9 - Gly 765
Purity	90%
Source	Rabbit LDL R, Fc Tag (orb668886) is expressed from human 293 cells (HEK293). It contains AA Ala 9 - Gly 765 (Accession # P20063-1).
MW	109.2 kDa
Expiration Date	6 months from date of receipt.