



Product Datasheet

PEA15 (phospho-Ser116) antibody (orb6689)

Description	PEA15 (phospho-Ser116) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PEA15 around the phosphorylation site of Ser116 (QP(p-S)EE)
Target	Phospho-PEA15 (Ser116)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15 kDa
Uniprot ID	Q15121
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...



Immunohistochemical staining of rat brain...