

Product Datasheet

SARS-CoV-2 Spike RBD Neutralizing Antibody (orb668900)

Description	SARS-CoV-2 Spike RBD Neutralizing Antibody
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Spike RBD
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Note	For research use only
Application notes	This antibody can be paired with other Anti-SARS-CoV-2 Spike S1 antibodies to detect SARS-CoV-2 Spike S1 protein in sandwich ELISA or LFA assay.
Isotype	Human IgG1 Human Kappa
Clonality	Recombinant
Purity	95%
Expiration Date	12 months from date of receipt.