

Product Datasheet

SARS-CoV-2 Spike NTD Antibody (orb668905)

Description	SARS-CoV-2 Spike NTD Antibody
Species/Host	Human
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Spike NTD
Preservatives	PBS, pH 7.4
Form/Appearance	Powder
Storage	-20°C
Note	For research use only
Isotype	Human IgG1 Human Kappa
Clonality	Recombinant
Purity	95%
Source	Anti-SARS-CoV-2 Spike NTD Antibody, Chimeric mAb, Human IgG1 (AM121) (orb668905) is a chimeric monoclonal antibody combining the constant domains of the human IgG1 molecule with mouse variable regions.
Expiration Date	12 months from date of receipt.