

Product Datasheet

Human FGFR1c protein (orb668918)

Description Recombinant human Fibroblast Growth Factor Receptor 1 protein

Reactivity Human

Tested Applications Functional Studies, IA, WB

Immunogen Unconjugated

Preservatives PBS 20% Glycerol

Form/Appearance Liquid

Storage -20° or -70°C

Tag 10xHis at the C-terminus

Note For research use only

Protein Sequence 22-375


Purity > 95%

Source HEK293 cells

MW 40.7

Uniprot ID **P11362**

Expiration Date 6 months from date of receipt.



Human
FGFR1c
protein