

Product Datasheet

Human FGFR2b protein (orb668919)

Description Recombinant human Fibroblast Growth Factor Receptor 2 protein

Reactivity Human

Tested Applications Functional Studies, IA, WB

Immunogen Unconjugated

Preservatives PBS 20% Glycerol

Form/Appearance Liquid

Storage -20° or -70°C

Tag 10xHis at the C-terminus

Note For research use only

Protein Sequence 22-378


Purity > 95%

Source HEK293 cells

MW 41.2

Uniprot ID **P21802**

Expiration Date 6 months from date of receipt.



Human
FGFR2b
protein