

Product Datasheet

Human FGFR3c protein (orb668922)

Description Recombinant human Fibroblast Growth Factor Receptor 3 protein

Reactivity Human

Tested Applications Functional Studies, IA, WB

Immunogen Unconjugated

Preservatives PBS 20% Glycerol

Form/Appearance Liquid

Storage -20° or -70°C

Tag 10xHis at the C-terminus

Note For research use only

Protein Sequence 23-375


Purity > 95%

Source HEK293 cells

MW 39.7

Uniprot ID **P22607**

Expiration Date 6 months from date of receipt.



Human
FGFR3c
protein