

Product Datasheet

Human PLA2R1 protein (orb668924)

Description Recombinant human Secretory Phospholipase A2 Receptor 1 protein

Reactivity Human

Tested Applications Functional Studies, IA, WB

Immunogen Unconjugated

Preservatives PBS 20% Glycerol

Form/Appearance Liquid

Storage -20° or -70°C

Tag 10xHis at the C-terminus

Note For research use only

Protein Sequence 21-1397


Purity > 95%

Source HEK293 cells

MW 160.7

Uniprot ID **Q13018**

Expiration Date 6 months from date of receipt.



Human
PLA2R1
protein