

Product Datasheet

Human PLA2R1 protein (orb668926)

Description	Recombinant human Secretory Phospholipase A2 Receptor 1 protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Immunogen	Unconjugated
Preservatives	PBS 20% Glycerol
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	10xHis at the C-terminus
Note	For research use only
Protein Sequence	21-515
Purity	> 95%
Source	HEK293 cells
MW	58.5
Uniprot ID	Q13018
Expiration Date	6 months from date of receipt.