

Product Datasheet

Human sVAP-1 protein (orb668930)

Description	Recombinant human Soluble Vascular Adhesion Protein-1 protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Immunogen	Unconjugated
Preservatives	PBS 20% Glycerol
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	10xHis at the N-terminus
Note	For research use only
Protein Sequence	29-763
Purity	> 95%
Source	HEK293 cells
MW	83
Uniprot ID	Q16853
Expiration Date	6 months from date of receipt.