


Product Datasheet

Human APP alpha protein (orb668932)

Description	Recombinant human Amyloid Precursor Protein, alpha form, APP alpha protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Immunogen	Unconjugated
Preservatives	PBS 20% Glycerol
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	10xHis at the N-terminus
Note	For research use only
Protein Sequence	18-612
Purity	> 95%
Source	HEK293 cells
MW	69.2
Uniprot ID	P05067
Expiration Date	6 months from date of receipt.



Human APP
alpha
protein