


Product Datasheet

Human APP beta protein (orb668933)

Description	Recombinant human Amyloid Precursor Protein, beta form, APP beta protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Immunogen	Unconjugated
Preservatives	PBS 20% Glycerol
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	10xHis at the N-terminus
Note	For research use only
Protein Sequence	18-596
Purity	> 95%
Source	HEK293 cells
MW	67.3
Uniprot ID	P05067
Expiration Date	6 months from date of receipt.



Human APP
beta protein