


Product Datasheet

Human ATX protein (orb668934)

Description	Recombinant human Autotaxin, ATX protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Immunogen	Unconjugated
Preservatives	PBS 20% Glycerol
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	10xHis at the C-terminus
Note	For research use only
Protein Sequence	49-863
Purity	> 95%
Source	HEK293 cells
MW	96
Uniprot ID	Q13822
Expiration Date	6 months from date of receipt.



Human
ATX
protein