

Product Datasheet

Human SARS-CoV-2 trimeric soluble full-length Spike protein (orb668935)

Description

Recombinant human SARS-CoV-2 trimeric

Reactivity

Human

Tested Applications

Functional Studies, IA, WB

Immunogen

Unconjugated

Preservatives

PBS 20% Glycerol

Form/Appearance

Liquid

Storage

-20° or -70°C

Tag

8xHis at the C-terminus

Note

For research use only

Protein Sequence

16-1213

Purity

> 95%

Source




HEK293 cells

MW

137

Uniprot ID
P0DTC2
Expiration Date

6 months from date of receipt.


Human
SARS-CoV-
2 trimeric
soluble
full-l...

Human
SARS-CoV-
2 trimeric
soluble
full-l...

Human
SARS-CoV-
2 trimeric
soluble
full-l...