


Product Datasheet

EZH2 antibody (orb66894)

Description	Rabbit polyclonal to EZH2.Enhancer of zeste homolog 2 (EZH2) belongs to the Polycomb-group (PcG) fam...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human EZH2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	85
Uniprot ID	Q15910
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...